

Product datasheet for RC210750L2V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

galectin 9 (LGALS9) (NM_002308) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: galectin 9 (LGALS9) (NM_002308) Human Tagged ORF Clone Lentiviral Particle

Symbol: galectin 9

Synonyms: HUAT; LGALS9A

Mammalian Cell

Selection:

None

Vector: pLenti-C-mGFP (PS100071)

Tag: mGFP

ACCN: NM_002308

ORF Size: 969 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC210750).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 002308.3</u>

 RefSeq Size:
 1739 bp

 RefSeq ORF:
 972 bp

 Locus ID:
 3965

 UniProt ID:
 000182

 Cytogenetics:
 17q11.2

Domains: Gal-bind lectin

MW: 35.9 kDa







Gene Summary:

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The protein encoded by this gene is an S-type lectin. It is overexpressed in Hodgkin's disease tissue and might participate in the interaction between the H&RS cells with their surrounding cells and might thus play a role in the pathogenesis of this disease and/or its associated immunodeficiency. Multiple alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]