

Product datasheet for **RC210754L4V**

TACC3 (NM_006342) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TACC3 (NM_006342) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TACC3
Synonyms:	ERIC-1; ERIC1; maskin; Tacc4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006342
ORF Size:	2514 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210754).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006342.1
RefSeq Size:	2788 bp
RefSeq ORF:	2517 bp
Locus ID:	10460
UniProt ID:	Q9Y6A5
Cytogenetics:	4p16.3
Domains:	TACC
MW:	90.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the transforming acidic coiled-coil protein family. The encoded protein is a motor spindle protein that may play a role in stabilization of the mitotic spindle. This protein may also play a role in growth and differentiation of certain cancer cells. [provided by RefSeq, Nov 2011]