

Product datasheet for **RC210755L2V**

DSCR3 (VPS26C) (NM_006052) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DSCR3 (VPS26C) (NM_006052) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VPS26C
Synonyms:	DCRA; DSCR3; DSCRA
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006052
ORF Size:	891 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210755).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006052.1
RefSeq Size:	3252 bp
RefSeq ORF:	894 bp
Locus ID:	10311
UniProt ID:	O14972
Cytogenetics:	21q22.13
Domains:	Vps26
MW:	33 kDa



[View online »](#)

Gene Summary:

The region of chromosome 21 between genes CBR and ERG (CBR-ERG region), which spans 2.5 Mb on 21q22.2, has been defined by analysis of patients with partial trisomy 21. It contributes significantly to the pathogenesis of many characteristics of Down syndrome, including morphological features, hypotonia, and cognitive disability. The DSCR3 (Down syndrome critical region gene 3) gene is found in this region and is predicted to contain eight exons. DSCR3 is expressed in most tissues examined. [provided by RefSeq, Jul 2008]