

## Product datasheet for **RC210760L2V**

### DNA Ligase I (LIG1) (NM\_000234) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DNA Ligase I (LIG1) (NM_000234) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNA Ligase I
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000234
ORF Size:	2757 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210760).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000234.1</a>
RefSeq Size:	3083 bp
RefSeq ORF:	2760 bp
Locus ID:	3978
UniProt ID:	<a href="#">P18858</a>
Cytogenetics:	19q13.33
Domains:	DNA_ligase
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair, DNA replication, Mismatch repair, Nucleotide excision repair



[View online »](#)

**MW:** 101.6 kDa

**Gene Summary:** This gene encodes a member of the ATP-dependent DNA ligase protein family. The encoded protein functions in DNA replication, recombination, and the base excision repair process. Mutations in this gene that lead to DNA ligase I deficiency result in immunodeficiency and increased sensitivity to DNA-damaging agents. Disruption of this gene may also be associated with a variety of cancers. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]