

Product datasheet for **RC210784A1V**

Human C13orf30 (FAM216B) (NM_182508) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C13orf30 (FAM216B) (NM_182508) AAV Particle
Tag: Myc-DDK
Symbol: C13orf30
Synonyms: C13orf30
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210784 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGACAAAAGCTGGAAAAGACAACAAAAGCTTTGGAATGTTCCACAACCTCCTTTTATTCGAGTTCTCT
CCTCCATCTATGACACTTCCTTACTAAAGGCCCTCAACCAAGGGCAACAGCGCTACTTTTACAGCATTAT
GAGGATTTACAACCTCCAGGCCCCAGTGGGAGGCCCTGCAGACCCGCTACATTACAGCCTTCAGCACCAA
CAGCTGCTTGGCTATATTACTCAACGGGAAGCCTGTCTTATGCTCTTGTACTTAGAGATTCAACCAAGA
GAGCCTCAGCCAAGGTAGCTCCTCAAAGAACCATTCCCCGAAAACCTCAGCCATGACAAGAAGATGTCC
ATCAGTACTACCTGTATCTGTGGTTCTACCTAGGGCCCAAAGTAAAAGGCGCCAAGTGCTCAGGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210784 protein sequence
Red=Cloning site Green=Tags(s)

MGQNWKRQQLWNVPQLPFIKRVPSIYDTSLLKALNQGQRYFYSIMRIYNSRPQWEALQTRYIHSLSHQ
QLLGITQREALSYALVLRDSTKRASAKVAPQRTIPKTSAMTRRCPVLPVSVLPRAQSKRRQVLRN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_182508



[View online »](#)

ORF Size:	417 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_182508.1, NP_872314.1</u>
RefSeq Size:	3318 bp
RefSeq ORF:	420 bp
Locus ID:	144809
UniProt ID:	<u>Q8N7L0</u>
Cytogenetics:	13q14.11
MW:	16.4 kDa