

## Product datasheet for **RC210814L1V**

### **RCE1 (NM\_001032279) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RCE1 (NM_001032279) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RCE1
Synonyms:	FACE2; RCE1A; RCE1B
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001032279
ORF Size:	987 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210814).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001032279.1</a>
RefSeq Size:	1485 bp
RefSeq ORF:	678 bp
Locus ID:	9986
UniProt ID:	<a href="#">Q9Y256</a>
Cytogenetics:	11q13.2
Protein Families:	Protease, Transmembrane
MW:	35.8 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an integral membrane protein which is classified as a member of the metalloproteinase family. This enzyme is thought to function in the maintenance and processing of CAAX-type prenylated proteins. [provided by RefSeq, Jul 2008]