

## Product datasheet for **RC210824L3V**

### WDPCP (NM\_001042692) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WDPCP (NM_001042692) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDPCP
Synonyms:	BBS15; C2orf86; CHDTHP; CPLANE5; FRITZ; FRTZ
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001042692
ORF Size:	1761 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210824).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001042692.1</a>
RefSeq Size:	2511 bp
RefSeq ORF:	1764 bp
Locus ID:	51057
UniProt ID:	<a href="#">O95876</a>
Cytogenetics:	2p15
MW:	66.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a cytoplasmic WD40 repeat protein. A similar gene in frogs encodes a planar cell polarity protein that plays a critical role in collective cell movement and ciliogenesis by mediating septin localization. Mutations in this gene are associated with Bardet-Biedl syndrome 15 and may also play a role in Meckel-Gruber syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014]