

Product datasheet for **RC210825L1V**

IFNA5 (NM_002169) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | IFNA5 (NM_002169) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | IFNA5 |
| Synonyms: | IFN-alpha-5; IFN-alphaG; INA5; INFA5; IeIF G |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_002169 |
| ORF Size: | 567 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210825). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_002169.1 |
| RefSeq Size: | 700 bp |
| RefSeq ORF: | 570 bp |
| Locus ID: | 3442 |
| UniProt ID: | P01569 |
| Cytogenetics: | 9p21.3 |
| Protein Families: | Druggable Genome, Secreted Protein |



[View online »](#)

| | |
|--------------------------|---|
| Protein Pathways: | Antigen processing and presentation, Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Cytosolic DNA-sensing pathway, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Regulation of autophagy, RIG-I-like receptor signaling pathway, Toll-like receptor signaling pathway |
| MW: | 21.9 kDa |
| Gene Summary: | Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. [UniProtKB/Swiss-Prot Function] |