

Product datasheet for RC210834L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ADAMTS5 (NM_007038) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: ADAMTS5 (NM_007038) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: ADAMTS5

Synonyms: ADAM-TS 5; ADAM-TS5; ADAM-TS 11; ADAMTS-5; ADAMTS-11; ADAMTS11; ADMP-2

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clo


Sequence:

The ORF insert of this clone is exactly the same as(RC210834).

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_007038

ORF Size: 2790 bp





ADAMTS5 (NM_007038) Human Tagged Lenti ORF Clone - RC210834L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 007038.2, NP 008969.1

 RefSeq Size:
 9663 bp

 RefSeq ORF:
 2793 bp

 Locus ID:
 11096

 UniProt ID:
 Q9UNA0

 Cytogenetics:
 21q21.3

Domains: tsp_1, Reprolysin, ACR

Protein Families: Druggable Genome, Protease, Secreted Protein

MW: 101.7 kDa

Gene Summary: This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with


thrombospondin motifs) protein family. Members of the family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The encoded preproprotein is proteolytically processed to generate the mature enzyme. This enzyme contains two C-terminal TS motifs and functions as an aggrecanase that cleaves

aggrecan, a major proteoglycan of cartilage, and may mediate cartilage destruction in


osteoarthritis. [provided by RefSeq, Feb 2016]



Product images:



Circular map for RC210834L3



Double digestion of RC210834L3 using Sgfl and Rsrll $\,$