

## Product datasheet for **RC210855**

### **DSCR10 (NM\_148676) Human Tagged ORF Clone**

#### Product data:

Product Type: Expression Plasmids  
Product Name: DSCR10 (NM\_148676) Human Tagged ORF Clone  
Tag: Myc-DDK  
Symbol: DSCR10  
Mammalian Cell Selection: Neomycin  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
ORF Nucleotide Sequence: >RC210855 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCAGATTGTGCAGGGGTTCCCTGCGGATGCACCCCTGTGTGCCCTGATGTGGACCTGCTCCTTCTCTGC  
TTCCTGGGCTCCAGACAGAAACACCGTACCCCTGCACCAGCCTCTGCCTCAGCTCTTCACAGTCAGCACA  
CCCTCCGCTTCTGTGTCAGGGTGTCTCTGCTGAGTCTGGTTATGGGATCCCTTTCTGTGCAGAGCCATGC  
TCCCGTGTCACTGTGTGCACCTCAAGCTGTCCCGTGTGTATGCCTGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210855 protein sequence  
Red=Cloning site Green=Tags(s)  
MQIVQGFADAPLCAIMWTCFLLPGLQTETPYPTSLCLSSSQSAHPPLPVRVFAESGYGIPFCAEPC  
SRVTVCHLQAVPVCMPV

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Chromatograms: [https://cdn.origene.com/chromatograms/mk6382\\_e05.zip](https://cdn.origene.com/chromatograms/mk6382_e05.zip)

Restriction Sites: SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_148676

**ORF Size:** 261 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_148676.2](#), [NP\\_683517.1](#)


**RefSeq Size:** 863 bp

**RefSeq ORF:** 263 bp


Locus ID: 259234  
 Cytogenetics: 21q22.13  
 MW: 9.3 kDa

**Product images:**


Circular map for RC210855



Western blot validation of overexpression lysate (Cat# [LY407748]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC210855 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified DSCR10 protein (Cat# [TP310855]). The protein was produced from HEK293T cells transfected with DSCR10 cDNA clone (Cat# RC210855) using MegaTran 2.0 (Cat# [TT210002]).