

Product datasheet for **RC210868A1V**

Human RNF122 (NM_024787) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RNF122 (NM_024787) AAV Particle
Tag: Myc-DDK
Symbol: RNF122
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210868 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGACCCATTCCAGTGGTGAACGGGTGTTTCTGTGGCCTGGGACTGGTTAGCACCAACAAGTCCTGCT
 CGATGCCACCCATCAGTTTCCAGGACCTCCGCTCAACATCTATATGGTCATCTTCGGCACAGGCATCTT
 TGTCTTCATGCTCAGCCTTATCTTCTGCTGCTATTTTATCAGCAAAGTGCAGAACAGGCACAGAGTGAG
 CGATACGATATAAGGAGGTGGTGCTTAAAGGTGATGCCAAGAAGTTACAATTATATGGGCAGACCTGCG
 CAGTCTGTCTGGAAGACTTCAAGGGGAAGGATGAGTTAGGCGTGCTCCCGTGCCAACACGCCTTTCACCG
 CAAGTGTCTGGTGAAATGGCTGGAAGTTCGCTGTGTCTGCCCCATGTGTAACAAGCCATTGCTAGTCCC
 TCAGAGGCCACGCAGAACATTGGGATTCTATTGGATGAGCTGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210868 protein sequence
 Red=Cloning site Green=Tags(s)

MHPFQWCNGCFGLGLVSTNKSCMPPIISFQDLPLNIYMVIFGTGIFVFMLSLIFCCYFISKLRNQAQSE
 RYGYKEVVLKGDAAKLQLYGQTCVCLDFKKGDELGVLPQHAFHRKCLVKWLEVRVCVPMCNKPIASP
 SEATQNIIGILLDELV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_024787


[View online »](#)

ORF Size:	465 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024787.2, NP_079063.2</u>
RefSeq Size:	1872 bp
RefSeq ORF:	468 bp
Locus ID:	79845
UniProt ID:	<u>Q9H9V4</u>
Cytogenetics:	8p12
MW:	17.5 kDa