

Product datasheet for **RC210868L4V**

RNF122 (NM_024787) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNF122 (NM_024787) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF122
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024787
ORF Size:	465 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210868).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024787.2 , NP_079063.2
RefSeq Size:	1872 bp
RefSeq ORF:	468 bp
Locus ID:	79845
UniProt ID:	Q9H9V4
Cytogenetics:	8p12
Protein Families:	Druggable Genome, Transmembrane
MW:	17.5 kDa



[View online »](#)

Gene Summary:

The encoded protein contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions. The encoded protein is localized to the endoplasmic reticulum and golgi apparatus, and may be associated with cell viability. [provided by RefSeq, Jul 2013]