

Product datasheet for **RC210880**

PMS2 (NM_000535) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	PMS2
Synonyms:	HNPCC4; MLH4; MMRCS4; PMS2CL; PMSL2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC210880 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

 ATGGAGCGAGCTGAGAGCTCGAGTACAGAACCTGCTAAGGCCATCAAACCTATTGATCGGAAGTCAGTCC
ATCAGATTTGCTCTGGGCAGGTGGTACTGAGTCTAAGCACTGCGGTAAAGGAGTTAGTAGAAAACAGTCT
GGATGCTGGTGCCACTAATATTGATCTAAAGCTTAAGGACTATGGAGTGGATCTTATTGAAGTTTCAGAC
AATGGATGTGGGGTAGAAGAAGAAAACCTCGAAGGCTTAAGTCTGAAACATCACACATCTAAGATTCAAG
AGTTTGCCGACCTAACTCAGGTTGAACTTTTGGCTTTGGGGGGAAGCTCTGAGCTCACTTTGTGCACT
GAGCGATGTACCATTTTACCTGCCACGCATCGGCGAAGGTTGGAAGTCTGAGTGTGTTGATCACAAT
GGGAAAATTATCCAGAAAACCCCTACCCCGCCCCAGAGGGACCACAGTCAGCGTGCAGCAGTTATTTT
CCACACTACCTGTGCCCATAAAGGAATTTCAAAGGAATATTAAGAAGGAGTATGCCAAAATGGTCCAGGT
CTTACATGCATACTGTATCATTTTCAAGGAGTCCGTGTAAGTTGCACCAATCAGCTTGACAAGGAAAA
CGACAGCCTGTGGTATGCACAGGTGGAAGCCCCAGCATAAAGGAAAAATATCGGCTCTGTGTTGGGCAGA
AGCAGTTGCAAAGCCTCATTCTTTTGTTCAGCTGCCCCCTAGTACTCCGTGTGTGAAGAGTACGGTTT
GAGCTGTTTCGGATGCTCTGCATAATCTTTTACATCTCAGGTTTCATTTCAATGCACGCATGGAGTT
GGAAGGAGTTCAACAGACAGACAGTCTTTTCTTTATCAACCGGCGGCCCTGTGACCCAGCAAAGGTCTGCA
GACTCGTGAATGAGGTCTACCACATGTATAATCGACACCAAGTATCCATTTGTTGTTCTTAACATTCTGT
TGATTGAGAAATGCGTTGATATCAATGTTACTCCAGATAAAAGGCAAAATTTTGCTACAAGAGGAAAAGCTT
TTGTTGGCAGTTTTAAAGACCTCTTTGATAGGAATGTTTGATAGTGTATGCAACAAGCTAAATGTGAGTC
AGCAGCCACTGCTGGATGTTGAAGTAACCTTAATAAAAAATGCATGCAGCGGATTTGGAAAAGCCCATGGT
AGAAAAGCAGGATCAATCCCTTTCATTAAGGACTGGAGAAGAAAAAAGACGTGTCCATTTCCAGACTG
CGAGAGGCCTTTTCTCTCGTCACACAAGAGAAACAAGCCTCACAGCCCAAAGACTCCAGAACCAGAA
GGAGCCCTCTAGGACAGAAAAGGGGTATGCTGTCTTCTAGCACTTCAGGTGCCATCTCTGACAAAGGCGT
CCTGAGACCTCAGAAAGAGGCAGTGAGTCCAGTCACGGACCCAGTGACCCTACGGACAGAGCGGAGGTG
GAGAAGGACTCGGGGCACGGCAGCACTTCCGTGGATTCTGAGGGGTTGAGCATCCAGACACGGGCAGTC
ACTGCAGCAGCGAGTATGCGGCCAGTCCCCAGGGGACAGGGGCTCGCAGGAACATGTGGACTCTCAGGA
GAAAGCGCCTGAACTGACGACTCTTTTTCAGATGTGGACTGCCATTCAAACAGGAAGATACCGGATGT
AAATTTTCAGTTTTGCCTCAGCCAACTAATCTCGCAACCCCAAACACAAAGCGTTTTAAAAAAGAGAAA
TTCTTTCCAGTTCTGACATTTGTCAAAAGTTAGTAAATACTCAGGACATGTGAGCCTCTCAGGTTGATGT
AGCTGTGAAAATTAATAAGAAAGTTGTGCCCTGGACTTTTCTATAAGTTCTTTAGCTAAACGAATAAAG
CAGTTACATCATGAAGCACAGCAAAGTGAAGGGGAACAGAATTACAGGAAGTTTAGGGCAAAGATTTGTC
CTGGAGAAAATCAAGCAGCCGAAGATGAATAAGAAAAGAGATAAGTAAACGATGTTTGCAGAAATGGA
AATCATTTGGTCAGTTTAACTTGGGATTTATAATAACCAAAGTGAATGAGGATATCTTCATAGTGGACCA
CATGCCACGGACGAGAAGTATAACTTCGAGATGCTGCAGCAGCACACCGTGCTCCAGGGGCAGAGGCTCA
TAGCACCTCAGACTCTCAACTTAACTGCTGTTAATGAAGCTGTTCTGATAGAAAATCTGGAATATTTAG
AAAGAATGGCTTTGATTTTGTATCGATGAAAATGCTCCAGTCACTGAAAGGGCTAACTGATTTCTTTG
CCAAGTAGTAAAACTGGACCTTCGGACCCAGGACGTCGATGAACTGATCTTCATGCTGAGCGACAGCC
CTGGGGTCATGTGCCGGCCTTCCGAGTCAAGCAGATGTTTGCCTCCAGAGCCTGCCGGAAGTCGGTGTAT
GATTGGGACTGCTCTTAACACAAGCGAGATGAAGAACTGATCACCACATGGGGGAGATGGACCACCC
TGGAAGTGTCCCATGGAAGGCCAACCATGAGACACATCGCCAACCTGGGTGTCATTTCTCAGAAC

 ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210880 protein sequence
 Red=Cloning site Green=Tags(s)


MERAESSSTEPAKAIKPIDRKSVHQICSGQVVLSTAVKELVENSLDAGATNIDLKLDYGVDLIEVSD
 NGCGVEEENFEGTLKHHTSKIQEFADLTQVETFGFRGEALSSLCALSDVTISTCHASAKVGTRLMFDHN
 GKIIQKTPYPRRGTTVSQQLFSTLPVRHKEFQRNIKEYAKMVQLHAYCIIISAGIRVSCTNQLGQGK
 RQPVVCTGGSPSIKENIGSVFGQKQLQSLIPFVQLPPSDSVCEEYGLSCSDALHNLFYISGFISQCTHGV
 GRSSTRQFFFINRRPCDPAKVCRLVNEVYHMYNRHQYPFVVLNISVDSECVDINVTPDKRQILLQEEKL
 LLAVLKTSLIGMFDSDVNKLNVSQQPLLDVEGNLIKMHAADEKPMVEKQDQSPSLRTGEEKDVSISRL
 REAFSLRHTTENKPHSPKTPEPRRSPLGQKRGMLSSSTSGAISDKGVLRPQKEAVSSSHGPSDPTDRAEV
 EKDSGHGSTSVDSSEGFSDPTGSHCSSEYAASSPGDRGSQEHVDSQEKAPETDDSFSDVDCHSNQEDTGC
 KFRVLPQPTNLATPNTKRFKKEEILSSSDICQKLVNTQDMSASQVDVAVKINKKVPLDFSISSLAKRIK
 QLHHEAQSEGEQNYRKFRKICPGENQAAEDELKEISKTMFAEMEIIIGQFNLGFIITKLNEDIFIVDQ
 HATDEKYNFEMLQQHTVLQGGRLIAPQTLNLTAVNEAVLIENLEIFRKNGFDFVIDENAPVTERAKLISL
 PTSKNWTFGPQDVDELIFMLSDSPGVMCRPSRVKQMFASRACRKSVMIGTALNTSEMKKLITHMGEMDHP
 WNCPHGRPTMRHIANLGVISQN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6692_d12.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_000535


ORF Size: 2586 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.


The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_000535.7</u>
RefSeq Size:	2851 bp
RefSeq ORF:	2589 bp
Locus ID:	5395
UniProt ID:	<u>P54278</u>
Cytogenetics:	7p22.1
Domains:	DNA_mis_repair, HATPase_c
Protein Families:	Druggable Genome
Protein Pathways:	Mismatch repair
MW:	95.8 kDa
Gene Summary:	The protein encoded by this gene is a key component of the mismatch repair system that functions to correct DNA mismatches and small insertions and deletions that can occur during DNA replication and homologous recombination. This protein forms heterodimers with the gene product of the mutL homolog 1 (MLH1) gene to form the MutL-alpha heterodimer. The MutL-alpha heterodimer possesses an endonucleolytic activity that is activated following recognition of mismatches and insertion/deletion loops by the MutS-alpha and MutS-beta heterodimers, and is necessary for removal of the mismatched DNA. There is a DQHA(X)2E(X)4E motif found at the C-terminus of the protein encoded by this gene that forms part of the active site of the nuclease. Mutations in this gene have been associated with hereditary nonpolyposis colorectal cancer (HNPCC; also known as Lynch syndrome) and Turcot syndrome. [provided by RefSeq, Apr 2016]


Product images:




Circular map for RC210880



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PMS2 (Cat# RC210880, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PMS2 (Cat# [TA506658]). Positive lysates [LY424658] (100ug) and [LC424658] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY424658]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC210880 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PMS2 protein (Cat# [TP310880]). The protein was produced from HEK293T cells transfected with PMS2 cDNA clone (Cat# RC210880) using MegaTran 2.0 (Cat# [TT210002]).