

Product datasheet for RC210885

KRTAP13-2 (NM_181621) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: KRTAP13-2 (NM_181621) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: KRTAP13-2
Synonyms: KAP13-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC210885 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCTACAACCTGCTGCTGGAACTTCTCCTCCCGCTCCTGTGGTACTACCTGCGCTACCCAGCAT
 CCTCACGTGGCTTTTCTACCCAGCAATCTGGTCTACAGCACTGACCTCTGCTCTCCAGCACCTGCCA
 GCTGGGTCTCTCTATAGGGCTGTCAGGAGATCTGCTGGGAGCCCACCAGCTGCCAGAGCTCTAT
 GTGGAGTCCAGCCCTGCCAGACCTCTGCTACCGCCCCAGAACCTCCTGCTCTGCAGTCTTGAAGA
 CGACTTACTCTGGTCTCTAGGCTTTGGATCCAGCAGCTGCCGCTCCCTGGGCTATGGATCGAGGAGCTG
 CTACTCAGTGGGCTGTGGCTCCAGTGGTGTGAGATCCCTGGGTTATGGAAGCTGTGGCTTCCCTTCCCTC
 GGCTATGGATCTGGATTCTGCCGCCAACCTACTTGGCTTCTAGGAGCTGCCAGTCTCCTGTTACAGAC
 CAGCCTATGGATCAACCTTCTGCAGATCAACTTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210885 protein sequence
 Red=Cloning site Green=Tags(s)

MSYNCCSGNFSSRSCGDYLRYPASSRGSYPSNLVYSTDLCSPSTCQLGSSLYRGCQEICWEPTSCQTSY
 VESSPCQTSCYRPRTSLLCSPCKTTYSGSLGFGSSCSRSLGYGSRSCYSVGCSSGVRSLGYGSCGFPSL
 YGSGFCRPTYLASRSCQSPCYRPAYGSTFCRSTC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6383_e05.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_181621

ORF Size: 525 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181621.4](#)

RefSeq Size: 849 bp

RefSeq ORF: 528 bp

Locus ID: 337959


UniProt ID: [Q52LG2](#)

Cytogenetics: 21q22.11


MW: 18.7 kDa

Gene Summary: In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC210885



Western blot validation of overexpression lysate (Cat# [LY405729]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC210885 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).