

Product datasheet for **RC210885L4V**

KRTAP13 (KRTAP13-2) (NM_181621) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | KRTAP13 (KRTAP13-2) (NM_181621) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | KRTAP13-2 |
| Synonyms: | KAP13-2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_181621 |
| ORF Size: | 525 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210885). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_181621.2 |
| RefSeq Size: | 849 bp |
| RefSeq ORF: | 528 bp |
| Locus ID: | 337959 |
| UniProt ID: | Q52LG2 |
| Cytogenetics: | 21q22.11 |
| MW: | 18.7 kDa |


[View online »](#)

Gene Summary:

In the hair cortex, hair keratin intermediate filaments are embedded in an interfilamentous matrix, consisting of hair keratin-associated proteins (KRTAP), which are essential for the formation of a rigid and resistant hair shaft through their extensive disulfide bond cross-linking with abundant cysteine residues of hair keratins. The matrix proteins include the high-sulfur and high-glycine-tyrosine keratins (By similarity).[UniProtKB/Swiss-Prot Function]