

## Product datasheet for **RC210888L3V**

### ING1 (NM\_198218) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ING1 (NM_198218) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ING1
Synonyms:	p24ING1c; p33; p33ING1; p33ING1b; p47; p47ING1a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198218
ORF Size:	840 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210888).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198218.1</a>
RefSeq Size:	2188 bp
RefSeq ORF:	633 bp
Locus ID:	3621
UniProt ID:	<a href="#">Q9UK53</a>
Cytogenetics:	13q34
Protein Families:	Druggable Genome, Transcription Factors
MW:	31.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a tumor suppressor protein that can induce cell growth arrest and apoptosis. The encoded protein is a nuclear protein that physically interacts with the tumor suppressor protein TP53 and is a component of the p53 signaling pathway. Reduced expression and rearrangement of this gene have been detected in various cancers. Multiple alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]