

## Product datasheet for **RC210891L4V**

### **S100A5 (NM\_002962) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	S100A5 (NM_002962) Human Tagged ORF Clone Lentiviral Particle
Symbol:	S100A5
Synonyms:	S100D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002962
ORF Size:	333 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210891).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002962.1</a>
RefSeq Size:	710 bp
RefSeq ORF:	279 bp
Locus ID:	6276
UniProt ID:	<a href="#">P33763</a>
Cytogenetics:	1q21.3
MW:	12.8 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein has a  $\text{Ca}^{2+}$  affinity 20- to 100-fold higher than the other S100 proteins studied under identical conditions. This protein also binds  $\text{Zn}^{2+}$  and  $\text{Cu}^{2+}$ , and  $\text{Cu}^{2+}$  strongly which impairs the binding of  $\text{Ca}^{2+}$ . This protein is expressed in very restricted regions of the adult brain. [provided by RefSeq, Jul 2008]