

Product datasheet for **RC210894A1V**

Human ETFB (NM_001985) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ETFB (NM_001985) AAV Particle
Tag: Myc-DDK
Symbol: ETFB
Synonyms: FP585; MADD
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210894 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGAGCTGCGCGTGCTCGTAGCTGTCAAGAGGGTCATCGACTACGCCGTGAAGATCCGAGTGAAGC
CTGACAGGACCGGTGTGGTCACGGATGGTGTGAAGCACTCCATGAACCCCTTCTGTGAGATCGCGGTGGA
GGAGGCTGTGCGGCTCAAGGAGAAGAAGCTGGTGAAGGAGGTATCGCCGTCAGCTGTGGCCTGCACAG
TGCCAGGAGACGATTCTGACGCCCTGGCCATGGGTGCAGACCGAGGTATCCACGTGGAGGTGCCCCAG
CAGAAGCAGAACGCTTGGTCCCCTGCAGGTGGCTCGGGTCTGGCCAAGCTGGCAGAGAAGGAGAAGGT
GGACCTGGTGCTGCTGGGCAAACAGGCCATCGATGATGACTGTAAACCAGACAGGGCAGATGACAGCTGGA
TTTCTTGACTGGCCACAGGGCACATTCCCTCCAGGTGACGCTGGAGGGGGACAAGTTGAAAGTGGAGC
GGGAGATCGATGGGGGCTGGAGACCCTGCGCCTGAAGCTGCCAGCTGTGGTGACAGCTGACCTGAGGCT
CAACGAGCCCCGCTACGCCACGCTGCCAACATCATGAAAGCCAAGAAGAAGATCGAGGTGATCAAG
CCTGGGACCTGGGTGTGGACCTGACCTCCAAGCTCTCTGTGATCAGTGTGGAGGACCCGCCACGCGCA
CGGCCGGCGTCAAGGTGGAGACCACTGAGGACCTGGTGCCAAGCTGAAGGAGATTGGGCGGATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210894 protein sequence Red=Cloning site Green=Tags(s) MAELRVLVAVKRVIDYAVKIRVKPDRTGVVTDGVKHSMPFCEIAVEEAVRLKEKKLVKEVIAVSCGPAQ CQETIRTALAMGADRGIVHVEVPPAEERLGPLQVARVLAKLAEKEKVDLVLLGKQAIDDDCNQTGQMTAG FLDWPQGTFFASQVTLEGDKLKVEREIDGGLETLRLLKLPVVTDLRLNEPRYATLPNIMKAKKKKIEVIK PGDLGVDLTSKLSVISVEDPPQRTAGVKVETTEDLVAKLKEIGRI TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001985
ORF Size:	765 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001985.2</u>
RefSeq Size:	933 bp
RefSeq ORF:	768 bp
Locus ID:	2109
UniProt ID:	<u>P38117</u>
Cytogenetics:	19q13.41
MW:	27.8 kDa