

Product datasheet for **RC210901L1V**

TNFSF18 (NM_005092) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | TNFSF18 (NM_005092) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | TNFSF18 |
| Synonyms: | AITRL; GITRL; hGITRL; TL6; TNLG2A |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_005092 |
| ORF Size: | 531 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210901). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_005092.3 |
| RefSeq Size: | 748 bp |
| RefSeq ORF: | 534 bp |
| Locus ID: | 8995 |
| UniProt ID: | Q9UNG2 |
| Cytogenetics: | 1q25.1 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Cytokine-cytokine receptor interaction |



[View online »](#)

MW: 20.31 kDa

Gene Summary: The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptor TNFRSF18/AITR/GITR. It has been shown to modulate T lymphocyte survival in peripheral tissues. This cytokine is also found to be expressed in endothelial cells, and is thought to be important for interaction between T lymphocytes and endothelial cells. [provided by RefSeq, Jul 2008]