

Human Tuberoinfundibular peptide (PTH2) (NM 178449) AAV Particle

Product Type:	AAV Particles
Product Name:	Human Tuberoinfundibular peptide (PTH2) (NM_178449) AAV Particle
Tag:	Myc-DDK
Symbol:	Tuberoinfundibular peptide
Synonyms:	TIP39
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210906 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGAGATCTGCCGCCGCGATCGCC ATGGAGACCCGCCAGGTGTCCAGGAGCCCTCGGGTTCGGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGGTGGTGCCTGGGGCGTCCGCACTGCCTCGGGAGTCGCCCTGCCCGCGGTCTGAGCTCCGCCCCAGGACGGGCTGGGCGGATCCCGCCACCCCAAGCCGCGGAGGAGCCTGGCGCTGGCGGACGACGCGGCC TTCCGGGAGCGCGCGGTTTGCTGGCCGCCCTCGAGCGCCGCCACTGGCTGAACCTCGTACATGCACAAGCTGCTGGTGTGGATGCGCCC ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC210906 protein sequence Red=Cloning site Green=Tags(s) METRQVSRSPRVRLLLLLLLLLLVVPWGVRTASGVALLPPVGLSLRPPGRAWADPATPRPRSLALADDAAFRERARLLAALERRHWLNSYMHKLLVLDAF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_178449
ORF Size:	300 bp

[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_178449.2
RefSeq Size:	459 bp
RefSeq ORF:	303 bp
Locus ID:	113091
UniProt ID:	Q96A98
Cytogenetics:	19q13.33
MW:	11.2 kDa