

## Product datasheet for **RC210929L2V**

### PHACTR1 (NM\_030948) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PHACTR1 (NM_030948) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHACTR1
Synonyms:	DEE70; dj257A7.2; EIEE70; RPEL; RPEL1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_030948
ORF Size:	1740 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210929).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_030948.1</a>
RefSeq Size:	2896 bp
RefSeq ORF:	1743 bp
Locus ID:	221692
UniProt ID:	<a href="#">Q9C0D0</a>
Cytogenetics:	6p24.1
MW:	66.3 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the phosphatase and actin regulator family of proteins. This family member can bind actin and regulate the reorganization of the actin cytoskeleton. It plays a role in tubule formation and in endothelial cell survival. Polymorphisms in this gene are associated with susceptibility to myocardial infarction, coronary artery disease and cervical artery dissection. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Apr 2016]