

## Product datasheet for **RC210941L4V**

### SDR O (SDR9C7) (NM\_148897) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SDR O (SDR9C7) (NM_148897) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SDR O
Synonyms:	ARCI13; RDHS; SDR-O; SDRO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_148897
ORF Size:	939 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210941).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_148897.1</a>
RefSeq Size:	1765 bp
RefSeq ORF:	942 bp
Locus ID:	121214
UniProt ID:	<a href="#">Q8NEX9</a>
Cytogenetics:	12q13.3
Protein Families:	Druggable Genome
MW:	35.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein with similarity to the short-chain dehydrogenase/reductase (SDR) family but has not been shown to have retinoid or dehydrogenase activities. [provided by RefSeq, Apr 2010]