

ORF Size:	405 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005478.3</u>
RefSeq Size:	718 bp
RefSeq ORF:	408 bp
Locus ID:	10022
UniProt ID:	<u>Q9Y5Q6</u>
Cytogenetics:	1p31.3
MW:	15.3 kDa