

## Product datasheet for **RC210973L2V**

### ADCK5 (NM\_174922) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ADCK5 (NM_174922) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADCK5
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_174922
ORF Size:	1740 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210973).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_174922.3</a>
RefSeq Size:	1969 bp
RefSeq ORF:	1743 bp
Locus ID:	203054
UniProt ID:	<a href="#">Q3MIX3</a>
Cytogenetics:	8q24.3
Protein Families:	Druggable Genome, Protein Kinase
MW:	65.7 kDa



[View online »](#)

**Gene Summary:**

The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).[UniProtKB/Swiss-Prot Function]