

## Product datasheet for **RC210976L3V**

### SLC9A3 (NM\_004174) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC9A3 (NM_004174) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC9A3
Synonyms:	DIAR8; NHE-3; NHE3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004174
ORF Size:	2502 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210976).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004174.1</a>
RefSeq Size:	2584 bp
RefSeq ORF:	2505 bp
Locus ID:	6550
UniProt ID:	<a href="#">P48764</a>
Cytogenetics:	5p15.33
Protein Families:	Druggable Genome, Transmembrane
MW:	92.7 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is an epithelial brush border Na/H exchanger that uses an inward sodium ion gradient to expel acids from the cell. Defects in this gene are a cause of congenital secretory sodium diarrhea. Pseudogenes of this gene exist on chromosomes 10 and 22. [provided by RefSeq, Mar 2016]