

Product datasheet for **RC211028**

PCDHGB2 (NM_032096) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHGB2 (NM_032096) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PCDHGB2
Synonyms:	PCDH-GAMMA-B2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC211028 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

 ATGAAAGCGAGCTCAGGGAGGTGCGGGCTGGTGCGGTGGCTGCAGGTACTGTTGCCCTTCCTGTTGTCTT
TGTTCCCGGGGCTCTCCAGTCCAGATCCGCTATTCAATTCCAGAGGAGCTGGCCAAAACTCGGTCGT
AGGAAACCTCGCCAAGGATCTGGGGCTCAGCGTCCGGGACTTGCCAGCCCGGAAGCTGCGGGTTAGCGCG
GAGAAGGAATATTTACAGTAACCCAGAAAGCGGAGACTTACTGTGAGTGACAGAATAGACCGAGAAC
AGATATGCGGGAAGCAGCCTCTGTGTCTGGATTTGATACTGTGCTGAAAACTCACTAAATATTTT
CTACATAGCAGTAATTGTGCAGGATATAATGATAATACCCGCTATTCAAACAGACTAAGATTAATTTA
AAAATTGGCGAATCCACTAAGCCAGGTACAACATTTCACTTGACCCAGCCCTGGATTAGATGTTGGTC
CTAACTCACTACAAAGATACCACCTTAATGACAACGAGTACTTTGATCTCGCTGAGAAACAGACTCCAGA
TGGTCGTAATATCCTGAGTTGATTCTAAACACTCTCTGGACAGAGAAGAGCACAGTTTACATCAATTG
GTCTCTACAGCTGTGGATGGCGGAGACCCACCTCAAAGTGGCACGACCCAAATCCGAATCAAAGTACCG
ATGCCAACGATAACCTCCAGTGTTCAGCCAGGACGTGTACAGGGTACACCTGAGGGAGGACGTGCCGCC
GGGCTTCTTTGTGCTTCAAGTGACAGCCACCGACCGGGATGAAGGCATAAACGCAGAGATCACCTACTCC
TTTCATAATGTGGACGAACAAGTAAACACTTTTTCACTTAAATGAAAAACAGGAGAAATCACGACAA
AGGATGATTTGGATTTTGAGATTGCAAGTAGTTACACTCTGAGTATCGAAGCAAAAGATCCTGGAGATCT
AGCAGCCCACTGCAGTATCCAAGTTGAAATCTTGATGACAACGATTGTGCACCTGAAGTTATTGTGACT
TCAGTATCTACTCCCCTACCGGAGGATTCGCCACCAGGAACAGTGATCGCCTTGATAAAACAGAGAGACA
GAGACTCTGGAGAAAATGGAGAAGTTTACTGCCAAGTGTGGGAAATGCCAAGTTATTTTGAAATCTTC
CTCAAAGAACTATTACAACTAGTGACAGACGGCGCTCTGGACCGGGAGGAGATCCAGAATACAATCTC
ACCATCACAGCCACCGACGGGGCAAGCCGCCCTCTCCTCCAGCATAATTGTACCCCTGCACATCTCCG
ACGTCAACGATAATGCCCGAGTTTCCAACAGACTTCTACATGGTTACGTGGCAGAGAACAATCCTCC
TGGCGCTCTATCGCTCAAATCAGTGCCTCTGACCCTGACTTGGGCCCAAGTGGCCAAAGTTCTACTCC
ATCGTAGCGAGCGACCTGAAGCCGCGGGAGATTTATCCTACGTGTCCGTGAGCGCGCAGAGCGGGTGG
TGTTTCGCGCAGCGCGCTTCGATCATGAGCAGCTGCGCGCTTCGAGCTCACACTGCAGGCCCGCGACCA
GGGCTCGCCCGCTCAGCGCAACGTGAGCCTGCGCGTGTAGTGGCGACCTCAATGACAATGCGCCA
CGGGTGTGTACCCCGCTGGGGCTGATGGCTCCGCCCTCTCGATATGGTGCCACGCGCCGAGAGC
CCGGCTACCTGGTGACCAAGGTGGTGGCGGTGGACGCACTCAGGACACAACGCTTGGCTGTCTACCA
CGTGCTGCAGGCCAGCGAGCCCGGGCTCTCAGCTGGGGTTGCGCACGGGTGAGGTGCGCACAGCGCT
GCCTTGGGCGACAGGGACGCGGCCCGCCAGCGCTGCTGGTGGTGTGCGTGATGGAGGACAGCCGCCAC
TCTCCGTACGGCCACGCTGCACCTAATCTTCGCGGATAGCCTGCAAGAGGTATTGCCAGACCTCAGCGA
CCGCCGGAGCCCTTGACCCCGAGGCAAACTGCAGTTTACCTGGTTGTGGCCTTGGCCTTGATCTCA
GTGCTCTTCTTCGCGGTGATTCTGGCAATCTCCCTGCGCCTGCGACTCTCTCCAGGTCAGATGCTT
GGGACTGTTTTAGCCTGGTCTCAGCTCCAAGCCTGGACCTGGGGTTCTCCCAATTACAGTGAGGGTAC
ATTGCCCTATTCTACAACCTGTGTGTTGCCTCACAATCAGCAAGACAGAGTTCAATTTCTGAACATA
ACCCCGGAATTGGTTCCCGCGCAAGATCTCGTCTGTGACAATGCCCTTTGGGAACAAAATACAAATCATG
GAGCCGCTGGGGTCCCTTTGCCTCAGATACTATTTGAAGGTGAGCTTTAAT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC211028 protein sequence

Red=Cloning site Green=Tags(s)

MKASSGRGCLVRWLQVLLPFLLSLFPGALPVQIRYSIPEELAKNSVVGNLAKDLGLSVRDLPAKLRVSA
EKEYFTVNPESGDLLVSDRIDREQICGKQPLCVLDFDTVAENPLNIFYIAVIVQDINDNTPLFKQTKINL
KIGESTKPGTTFFPLDPALDSVGPNSLQRYHLNDNEYFDLAEKQTPDGRKYPELIKHSLEDREEHSLHQL
VLTAVDGGDPPQSGTTQIRIKVTDANDNPPVFSQDVYRVTLREDVPPGFFVLQVTATDRDEGINAEITYS
FHNVDEQVKHFFNLNEKTGEITTKDDLDFEIASSTLSIEAKDPGDAAHCSIQVEILDNDNCAPEVIVT
SVSTPLPEDSPPGTVIALIKTRDRDSGENGEVYCQVLGNAKFILKSSSKNYYKLVTGALDREEIPEYNL
TITATDGGKPPLSSSIIVTLHISDVNDNAPVFQQTSMVHVAENPPGASIAQISASDPDLGPSGQVSY
IVASDLKPREILSYVSVSAQSGVVFAQRAFDHEQLRAFELTLQARDQGSPALSANVSLRVLVGDLNDNAP
RVLYPALGPDGSALFDMVPRAAEPGYLVTKVAVDADSGHNAWLSYHVLQASEPGLFSLGLRTGEVARTAR
ALGDRDAARQLLVAVRDGGQPPLSATATLHLIFADSLQEVLPDLSDRREPSDPQAKLQFYLVALALIS
VLFFLAVILAIISLRRLSSRSDAWDCFQPLSSKPGPGVLPNYSEGTLPSYNLCVASQSAKTEFNFLNI
TPELVPAQDLVCDNASWEQNTNHAAGVPFASDTILKVSFN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk6692_f03.zip

Restriction Sites:

Sgfl-MluI


Cloning Scheme:




ACCN:	NM_032096
ORF Size:	2433 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_032096.1, NP_115267.1</u>
RefSeq Size:	2436 bp
RefSeq ORF:	2436 bp
Locus ID:	56103
UniProt ID:	<u>Q9Y5G2</u>
Cytogenetics:	5q31.3
Domains:	CA
MW:	88.4 kDa
Gene Summary:	This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]


Product images:



Circular map for RC211028



Western blot validation of overexpression lysate (Cat# [LY410378]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC211028 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PCDHGB2 protein (Cat# [TP311028]). The protein was produced from HEK293T cells transfected with PCDHGB2 cDNA clone (Cat# RC211028) using MegaTran 2.0 (Cat# [TT210002]).