

Product datasheet for RC211037A1V

Human IL36 beta (IL36B) (NM_173178) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human IL36 beta (IL36B) (NM_173178) AAV Particle
Tag:	Myc-DDK
Symbol:	IL36 beta
Synonyms:	FIL1; FIL1-(ETA); FIL1H; FILI-(ETA); IL-1F8; IL-1H2; IL1-ETA; IL1F8; IL1H2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211037 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGAACCCACAACGGGAGGCAGCACCCAAATCCTATGCTATTCGTGATTCTCGACAGATGGTGTGGGTCC TGAGTGGAAATCTTTAATAGCAGCTCCTCTTAGCCGCAGCATTAAAGCCTGTCACTCTTCATTTAATAGC CTGTAGAGACACAGAATTCAGTGACAAGGAAAGGGTAATATGGTTTACCTGGGAATCAAGGGAAAAGAT CTCTGTCTCTTCTGTGCAGAAATTCAGGGCAAGCCTACTTTCAGCTTAAGGAAAAAATATCATGGACC TGTATGTGGAGAAGAAAGCACAGAAGCCCTTCTCTTTTCCACAATAAAGAAGGCTCCACTTCTGTCTT TCAGTCAGTCTCTTACCCTGGCTGGTTCATAGCCACCTCCACCACATCAGGACAGCCCATCTTCTCACC AAGGAGAGAGGCATAACTAATAACACTAACTTCTACTTAGATTCTGTGGAA ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC211037 protein sequence Red=Cloning site Green=Tags(s) MNPQREAPKSYAIRDSRQMVVLSGNSLIAAPLSRSIKPVTLHLIACRDTEFSDKEKGNMYYLGIKGD LCLFCAEIQGKPTLQLKEKNIMDLVVEKKAQKPFLFFHNKEGSTSVFQSVSPGWFIASTTSGQPIFLT KERGITNNTNFYLDSE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2


[View online »](#)

ACCN:	NM_173178
ORF Size:	471 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173178.1</u>
RefSeq Size:	1068 bp
RefSeq ORF:	474 bp
Locus ID:	27177
UniProt ID:	<u>Q9NZH7</u>
Cytogenetics:	2q14.1
MW:	17.7 kDa