

## Product datasheet for **RC211039L3V**

### IL17C (NM\_013278) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL17C (NM_013278) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IL17C
Synonyms:	CX2; IL-17C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013278
ORF Size:	591 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211039).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013278.3</a>
RefSeq Size:	1048 bp
RefSeq ORF:	594 bp
Locus ID:	27189
UniProt ID:	<a href="#">Q9P0M4</a>
Cytogenetics:	16q24.2
Protein Families:	Druggable Genome, Secreted Protein
MW:	21.8 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a T cell-derived cytokine that shares the sequence similarity with IL17. This cytokine was reported to stimulate the release of tumor necrosis factor alpha and interleukin 1 beta from a monocytic cell line. The expression of this cytokine was found to be restricted to activated T cells. [provided by RefSeq, Jul 2008]