

## Product datasheet for **RC211048L4V**

### **PAGE5 (NM\_001013435) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PAGE5 (NM_001013435) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAGE5
Synonyms:	CT16; CT16.1; CT16.2; GAGEE1; PAGE-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001013435
ORF Size:	330 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211048).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001013435.1</a>
RefSeq Size:	560 bp
RefSeq ORF:	333 bp
Locus ID:	90737
UniProt ID:	<a href="#">Q96GU1</a>
Cytogenetics:	Xp11.21
MW:	11.8 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of family of proteins that are expressed in a variety of tumors and in some fetal and reproductive tissues. The encoded protein may protect cells from programmed cell death. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found. [provided by RefSeq, Jan 2015]