

Product datasheet for **RC211076L3V**

UBL5 (NM_001048241) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UBL5 (NM_001048241) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UBL5
Synonyms:	HUB1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001048241
ORF Size:	219 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211076).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001048241.1
RefSeq Size:	469 bp
RefSeq ORF:	222 bp
Locus ID:	59286
UniProt ID:	Q9BZL1
Cytogenetics:	19p13.2
MW:	8.5 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of a group of proteins similar to ubiquitin. The encoded protein is not thought to degrade proteins like ubiquitin but to affect their function through being bound to target proteins by an isopeptide bond. The gene product has been studied as a link to predisposition to obesity based on its expression in *Psammomys obesus*, the fat sand rat, which is an animal model for obesity studies. Variation in this gene was found to be significantly associated with some metabolic traits (PMID: 15331561) but not associated with childhood obesity (PMID: 19189687). Pseudogenes of this gene are located on chromosomes 3, 5 and 17. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2013]