

Product datasheet for **RC211086L4V**

RFX7 (NM_022841) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RFX7 (NM_022841) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RFX7
Synonyms:	RFXDC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022841
ORF Size:	4380 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211086).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022841.5 , NP_073752.5
RefSeq Size:	7203 bp
RefSeq ORF:	4383 bp
Locus ID:	64864
UniProt ID:	Q2KHR2
Cytogenetics:	15q21.3
Protein Families:	Transcription Factors
MW:	157.1 kDa



[View online »](#)

Gene Summary:

RFX7 is a member of the regulatory factor X (RFX) family of transcription factors (see RFX1, MIM 600006) (Aftab et al., 2008 [PubMed 18673564]).[supplied by OMIM, Mar 2009]