

Product datasheet for **RC211096A1V**

Human C11orf52 (NM_080659) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C11orf52 (NM_080659) AAV Particle
Tag: Myc-DDK
Symbol: C11orf52
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC211096 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGAAACCGGGTCTGCTGCGGAGGAAGCTGGAGCTGCCCATCAACTTTCCAGAAGAAAAAGAAAGAG
 GAAGCCAAACAAGACGGACACTGAAGCCGCAGCCACAACAGCTGCAGCAGAATCTCCAAAGGGCCATGA
 AACACAGGACATACGTATGAACGGGTGTACAGCAGCAAGGGTCTCAAGAGAGGAGTCCAGGCCTCATG
 TCGGAAGACAGCAACTTACATTATGCTGACATTCAAGTGTGCAGCCGTCCCCATGCCCGGAAGTGAAC
 ACGTGCATTTAGAAAACGCTACAGAGTATGCGACCCTTCGCTTCCCCAGGCCACACCTCGCTATGACAG
 CAAGAACGGGACCCTGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211096 protein sequence
 Red=Cloning site Green=Tags(s)

MGNRVCCGGSWSCPSTFQKKKKRGSQTRRTLKPQPQLQQNLPGHETTGHYERVLQQQSQERSPLM
 SEDSNLHYADIQCSRPHAREVKHVHLENATEYATLRFPQATPRYDSKNGTLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_080659
ORF Size: 369 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080659.1</u> , <u>NP_542390.1</u>
RefSeq Size:	1168 bp
RefSeq ORF:	372 bp
Locus ID:	91894
UniProt ID:	<u>Q96A22</u>
Cytogenetics:	11q23.1
MW:	14 kDa