

Product datasheet for **RC211100A1V**

Human NECAP1 (NM_015509) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NECAP1 (NM_015509) AAV Particle
Tag:	Myc-DDK
Symbol:	NECAP1
Synonyms:	DEE21; EIEE21
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211100 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGACCGAGTTGGAGTACGAGTCTGTGCTGTGTGAAGCCAGACGTCAGCGTCTACCGGATTCCGC
CCCGGGCCTCCAACCGCGGTTACAGGGCCTCTGACTGGAATTAGACCAGCCTGATTGGACTGGTCGCCT
CCGAATCACTTCAAAAGGAAGACTGCCTATATCAAACGAGGATAAAGTTTCAGGGGAGCTCTTTGCT
CAGGCACCACTAGAACAAATATCCTGGTATTGCTGTGGAGACGGTGACAGATTCTAGCCGCTACTTTGTAA
TCCGGATCCAGGATGGTACTGGGCGCAGTGCTTTCATTGGCATTGGCTTCACAGATCGGGGAGATGCCTT
CGACTTTAATGTCTCCTTGCAGGATCACTTCAAGTGGTAAAGCAGGAATCTGAGATTTCCAAGGAATCT
CAAGAAATGGATGCTCGTCCTAAGTTGGATCTGGGCTTCAAGGAAGGACAAACCATCAAGTTGTGTATCG
GGAACATTACAAACAAGAAAGGAGGTGCTTCTAAGCCCAGGACTGCAAGGGGTGGGGTCTGAGCTTACT
CCCACCCCGCCAGGAGGCAAGTCACTATTCCCCACCATCCTCCTCAGTTGCCATCAGCAATCATGTC
ACCCACCAACCATTCGAAATCTAACCATGGAGGCAGTGATGCAGATATCCTTTTAGATTGGATTCTC
CTGCTCCTGTACGACACCAGCACCAACTCCAGTTTCTGTAAGCAATGACTTGTGGGGAGACTTCAGCAC
TGCTCCAGCTCTGTTCCAAACAGGCACACAGCCATCCAACCTGGGTCCAGTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211100 protein sequence Red=Cloning site Green=Tags(s) MATELEYESVLCVKPDVSYYRIPPRASNRGYRASDWKLDQPDWTGRLRITSKGKTAYIKLEDKVSSELF QAPVEQYPGIAVETVTDSSRYFVIRIQDGTGRSAFIGIGFTDRGDAFDNFVSLQDHFQWVKQSEISKES QEMDARPKLDLGFKEGQTIKLCIGNITNKKGGASKPRTARGGGLSLPPPPGGKVTIPPPSSVAISNHV TPPPIPKSNHGGSDADILLDLSPAPVTPAPTPVSVSNDLWGDFSTASSVPNQAPQPSNVVQF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015509
ORF Size:	825 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015509.2</u>
RefSeq Size:	2600 bp
RefSeq ORF:	828 bp
Locus ID:	25977
UniProt ID:	<u>Q8NC96</u>
Cytogenetics:	12p13.31
MW:	29.7 kDa