

## Product datasheet for **RC211101L4V**

### GCET1 (SERPINA9) (NM\_175739) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GCET1 (SERPINA9) (NM_175739) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GCET1
Synonyms:	GCET1; SERPINA11; SERPINA11b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_175739
ORF Size:	1305 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_175739.3</a> , <a href="#">NP_783866.2</a>
RefSeq Size:	1851 bp
RefSeq ORF:	1254 bp
Locus ID:	327657
UniProt ID:	<a href="#">Q86WD7</a>
Cytogenetics:	14q32.13
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	48.6 kDa



[View online »](#)

**Gene Summary:**

Protease inhibitor that inhibits trypsin and trypsin-like serine proteases (in vitro). Inhibits plasmin and thrombin with lower efficiency (in vitro).[UniProtKB/Swiss-Prot Function]