

Product datasheet for **RC211103A1V**

Human CCDC28A (NM_015439) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CCDC28A (NM_015439) AAV Particle
Tag:	Myc-DDK
Symbol:	CCDC28A
Synonyms:	C6orf80; CCRL1AP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211103 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGAGGGCGGAGCCAAGAGCGACACTGGGGGAGCAGGAAAAGCGGGGCTTCCGCTTGGGGCATGGA
GGCTGTACCTCTTACGTCACTTCCGTAACAAACGGAGCTGCGGAGGAGCGGGTCCCGGGATGTGACCGG
GGCTCTGCTTGTGGCTGCGGCGGTGGCTTCTGAGGCTGTGGGTCTTTCGGGTGCGGAAGGGGGCCCC
AATACCTTCTTCTTCAGGTCTTAAGAAGCTGGCCGTGGTGAATAAGGAACCTAAAACAATGGAAGAGC
GGAAAGTGAAGAGGAGGAGTCCTAAGTCTTTAGTGCCCACTGTACTCAGTTGTCAATGCCAAAAAAA
TGCCATTCCAGTGAGTAAAAGCACAGGGTTTTCAAATCCTGCATCACAGTCAACTTCACAGCGACCAAAG
TTAAAAAGAGTGATGAAAGAAAAGACCAACCTCAGGGTGGAGAGGGCAAAGGCGCTCAGTCAACTCCGA
TCCAGCACTCCTTCTCACTGATGTCTCAGATGTTGAGGAGATGGAGAGAGGGCTGCTCAGTCTTTTGA
TGATTTCCACTCTGGAAGAACTTCAAGCATTTGGAATGAATGTTCCATTGAACAGATGGAACATGTTGCG
GGAATGCAGGAGAAATTAGCTCGCTTGAATTTGGAGCTCTATGGGAGTTAGAGGAACCTCCTGAGGATA
AGAGAAAAACAGCCAGTCACTCAATCTGGATAGGCTTCTGTGAGATTTAGAAGAATTGAATTCTCCAT
ACAAAACTCCATTTGGCAGATGCACAAGATGTTCCAAATACTTCTGCTAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211103 protein sequence Red=Cloning site Green=Tags(s) MPRAEPRATLGEQEKAGLPLGAWRLYLLRHFRKQTELRRSGSRDVTGALLVAAAVASEAVGSLRVAEGGP NTLLLQVLRSWPWCNKELKTMEERKVKRRSPKSFSAHCTQVVNAKKNAIPVSKSTGFSNPASQSTSQRPK LKRVMEKTKPQGEGKGAQSTPIQHSFLTVDSDVQEMERGLLSLLNDFHSGKLQAFGNECSIEQMEHVR GMQEKLARLNLELYGELEELPEDKRKTASDSNLDRLSDLEELNSSIQKLHLADAQDVPNTSAS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015439
ORF Size:	822 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015439.2, NP_056254.1</u>
RefSeq Size:	1504 bp
RefSeq ORF:	555 bp
Locus ID:	25901
UniProt ID:	<u>Q8IWP9</u>
Cytogenetics:	6q24.1
MW:	30.4 kDa