

Product datasheet for **RC211128L3V**

FAHD2A (NM_016044) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAHD2A (NM_016044) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAHD2A
Synonyms:	CGI-105
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016044
ORF Size:	942 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211128).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016044.1
RefSeq Size:	1313 bp
RefSeq ORF:	945 bp
Locus ID:	51011
UniProt ID:	Q96GK7
Cytogenetics:	2q11.1
Domains:	FAA_hydrolase
MW:	34.6 kDa



[View online »](#)

Gene Summary: May have hydrolase activity.[UniProtKB/Swiss-Prot Function]