

## Product datasheet for **RC211129L1V**

### IRX4 (NM\_016358) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IRX4 (NM_016358) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IRX4
Synonyms:	IRXA3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016358
ORF Size:	1557 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211129).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016358.1</a>
RefSeq Size:	2212 bp
RefSeq ORF:	1560 bp
Locus ID:	50805
UniProt ID:	<a href="#">P78413</a>
Cytogenetics:	5p15.33
MW:	54.3 kDa



[View online »](#)

**Gene Summary:** Likely to be an important mediator of ventricular differentiation during cardiac development.  
[UniProtKB/Swiss-Prot Function]