

## Product datasheet for **RC211134L1V**

### **INPP4B (NM\_003866) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	INPP4B (NM_003866) Human Tagged ORF Clone Lentiviral Particle
Symbol:	INPP4B
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003866
ORF Size:	2772 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211134).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003866.1</a>
RefSeq Size:	2937 bp
RefSeq ORF:	2775 bp
Locus ID:	8821
UniProt ID:	<a href="#">O15327</a>
Cytogenetics:	4q31.21
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system
MW:	104.6 kDa



[View online »](#)

**Gene Summary:**

INPP4B encodes the inositol polyphosphate 4-phosphatase type II, one of the enzymes involved in phosphatidylinositol signaling pathways. This enzyme removes the phosphate group at position 4 of the inositol ring from inositol 3,4-bisphosphate. There is limited data to suggest that the human type II enzyme is subject to alternative splicing, as has been established for the type I enzyme. [provided by RefSeq, Jul 2008]