

Product datasheet for **RC211152A1V**

Human **BCL2L2 (NM_004050) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Human BCL2L2 (NM_004050) AAV Particle
Tag: Myc-DDK
Symbol: BCL2L2
Synonyms: BCL-W; BCL2-L-2; BCLW; PPP1R51
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC211152 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGACCCAGCCTCGGCCCCAGACACACGGGCTCTGGTGGCAGACTTTGTAGGTTATAAGCTGAGGC
AGAAGGGTTATGTCTGTGGAGCTGGCCCCGGGAGGGCCAGCAGCTGACCCGCTGCACCAAGCCATGCG
GGCAGCTGGAGATGAGTTCGAGACCCGCTTCGGCGCACCTTCTCTGATCTGGCGGCTCAGCTGCATGTG
ACCCAGGCTCAGCCCAACACGCTTCACCCAGGTCTCCGATGAATTTTCAAGGGGGCCCAACTGGG
GCCGCTTGTAGCCTTCTTTGTCTTTGGGGCTGCACTGTGTGCTGAGAGTGTCAACAAGGAGATGGAACC
ACTGGTGGGACAAGTGCAGGAGTGGATGGTGGCTACCTGGAGACGCGGCTGGCTGACTGGATCCACAGC
AGTGGGGGCTGGGCGGAGTTCACAGCTCTATACGGGGACGGGGCCCTGGAGGAGGCGCGGCTCTGCGGG
AGGGGAAGTGGGCATCAGTGAGGACAGTGTGACGGGGGCCGTGGCACTGGGGGCCCTGGTAAGTGTAGG
GGCCTTTTGTAGCAAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211152 protein sequence
Red=Cloning site Green=Tags(s)

MATPASAPDTRALVADFVGYKLQKGYVCGAGPGGPAADPLHQAMRAAGDEFETRFRRTFSDLAAQLHV
TPGSAQQRFTQVSDLFQGGPNWGRLVAFFVFGAALCAESVNKEMEPLVGQVQEWVMVAYLETRLADWIHS
SGGWAEFTALYDGALEEARRLREGNWSVRTVLTGAVALGALVTVGAFASK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_004050
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004050.2</u>
RefSeq Size:	3621 bp
RefSeq ORF:	582 bp
Locus ID:	599
UniProt ID:	<u>Q92843</u>
Cytogenetics:	14q11.2
MW:	20.8 kDa