

Product datasheet for **RC211162A1V**

Human OR7A5 (NM_017506) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR7A5 (NM_017506) AAV Particle
Tag:	Myc-DDK
Symbol:	OR7A5
Synonyms:	HTPCR2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211162 representing NM_017506 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAACCAGGAAATGATACACAAATTCAGAATTTCTTCTTCTGGGATTTTCACAAGAACCTGGACTGC
AACCTTCCTCTTTGGGCTGTTCTGTCCATGTACCTGGTCACTGTGCTCGGGAACCTGCTCATCATCCT
GGCCACAATCTCAGACTCCACCTCCACACCCCATGTACTTCTTCTCTCAACCTGTCCTTTGCTGAC
ATTTGTGTTACTTCCACCACCATTCAAAAATGCTGATGAACATCCAGACACAGAACAAGTCATCACCT
ACATAGCCTGCCTCATGCAGATGTATTTTTCATACTCTTTGCTGGATTTGAAAACCTTCCTCTGTCCTG
GATGGCCTATGACCGGTTTGTGGCCATCTGTACCCCTGCCTACATGGTCATTATGAACCCCTCACCTC
TGTGGACTGCTGGTTCTAGCATCCTGGACCATGAGTGCTCTGTATTCTTGCTACAAATCTTAATGGTAG
TACGGCTGTCCTTCTGCACAGCCTTAGAAATCCCCCACTTTTCTGTGAACCTAATCAGGTCATCCAAC
TGCTTGTCTGATAGCTTTCTTAATCACATGGTGATATATTTTACAGTTGCGCTGCTGGGTGGAGGTCCC
CTGACTGGGATCCTTTACTCTTACTCTAAGATAATTTCTTCCATACATGCAATCTCATCAGCTCAGGGGA
AGTACAAGGCATTTTCCACCTGTGCATCTCACCTCTCAGTTGTCTCCTTATTTTATGGTGAATCCTAGG
GGGTACCTTAGTTCTGCTGCCACCCGCAACTCACACTCAAGTGCAACAGCCTCAGTGATGTACACTGTG
GTCACCCCATGCTGAACCCCTTATCTATAGTCTGAGGAATAAAGACATAAAGAGGGCTCTGGGAATAC
ATTTGTTGTGGGAACAATGAAAGGGCAATTTTCAAGAAGTGCCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211162 representing NM_017506 Red=Cloning site Green=Tags(s) MEPGNDTQISEFLLLGSQEPGLQPFLFGLFLSMYLVTVLGNLLIILATISDSHLHTPMYFFLSNLSFAD ICVTSTTIPKMLMNIQTQNKVITYIACLMQMYFFILFAGFENFLLSVMAYDRFVAICHPLHYMVMNPHL CGLLVLASWTMSALYSLLQILMVVRLSFCTALEIPHFFCELNQVIQLACSDSFLNHMVIYFTVALLGGGP LTGILYSYSKIISSIIHAISSAQGYKAFSTCASHLSVVSFLFYGAILGVYLSSAATRNSHSSATASVMYTV VTPMLNPFIIYSLRNKDIKRALGIHLLWGTMTKGQFFKKCP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_017506
ORF Size:	957 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017506.1</u>
RefSeq Size:	2138 bp
RefSeq ORF:	960 bp
Locus ID:	26659
UniProt ID:	<u>Q15622</u>
Cytogenetics:	19p13.12
MW:	35.4 kDa