

## Product datasheet for **RC211177L4V**

### ZIC1 (NM\_003412) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZIC1 (NM_003412) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZIC1
Synonyms:	BAIDCS; CRS6; ZIC; ZNF201
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003412
ORF Size:	1341 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003412.3</a>
RefSeq Size:	5231 bp
RefSeq ORF:	1344 bp
Locus ID:	7545
UniProt ID:	<a href="#">Q15915</a>
Cytogenetics:	3q24
Protein Families:	Druggable Genome
MW:	48.3 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. Members of this family are important during development. Aberrant expression of this gene is seen in medulloblastoma, a childhood brain tumor. This gene is closely linked to the gene encoding zinc finger protein of the cerebellum 4, a related family member on chromosome 3. This gene encodes a transcription factor that can bind and transactivate the apolipoprotein E gene. [provided by RefSeq, Jul 2008]