

Product datasheet for RC211187

Histone H1.3 (HIST1H1D) (NM_005320) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Histone H1.3 (HIST1H1D) (NM_005320) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Histone H1.3
Synonyms:	H1.3; H1D; H1F3; H1s-2; HIST1H1D
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC211187 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCGGAGACTGCTCCACTTGCTCCTACCATTCTGCACCCGCAGAAAAACACCTGTGAAGAAAAAGG
CGAAGAAGGCAGGCGCACTGCTGGGAAACGAAAGCATCCGGACCCCAAGTATCTGAGCTTATCACCAA
GGCAGTGGCAGCTTCTAAGGAGCGCAGCGCGTTTCTCTGGCCGCGCTTAAGAAAGCGTTGCGGCTGCC
GGCTACGATGTAGAAAAAACAACAGCCGATCAAGCTTGGCCTCAAGAGCTTGGTGAAGCAAGGTACCC
TGGTGCAGACAAAGGTACCGGTGCTTCTGGCTCCTTCAAACCAAGAAAGCGGCTTCCGGGAAGG
CAAACCAAGGCCAAAAAGGCTGGCGCAGCCAAGCCTAGGAAGCCTGCTGGGGCAGCCAAGAAGCCCAAG
AAGGTGGCTGGCGCCGCTACCCGAAGAAAAGCATCAAAAAGACTCCTAAGAAGGTAAAGAAGCCAGCAA
CCGCTGCTGGGACCAAGAAAGTGGCCAAGAGTGCAGAAAAGGTGAAAACACCTCAGCCAAAAAAGCTGC
CAAGAGTCCAGCTAAGGCCAAAGCCCCTAAGCCCAAGGCGGCCAAGCCTAAGTCGGGGAAGCCGAAGGT
ACAAAGGCAAAGAAGGCAGCTCCGAAGAAAAAG

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC211187 protein sequence
Red=Cloning site Green=Tags(s)

MSETAPLAPTIPAPAIEKTPVKKKAKKAGATAGKRKASGPPVSELITKAVAASKERSGVSLAALKKALAAA
 GYDVEKNNSRIKLGKSLVSKGTLVQTKGTGASGSFKLNKKAASGEGPKAKKAGAAKPRKPAGAAKPK
 KVAGAATPKKSIKKTTPKVKKPATAAGTKKVAKSARKVTPQPKAAKSPAKAKPKAAKPKSGKPKV
 TKAKKAAPKKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6373_g09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_005320

ORF Size: 663 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005320.3</u>
RefSeq Size:	777 bp
RefSeq ORF:	666 bp
Locus ID:	3007
UniProt ID:	<u>P16402</u>
Cytogenetics:	6p22.2
MW:	22.3 kDa
Gene Summary:	Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6. [provided by RefSeq, Aug 2015]

Product images:



Circular map for RC211187



Western blot validation of overexpression lysate (Cat# [LY417383]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC211187 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).