

## Product datasheet for **RC211202L4V**

### GLIS1 (NM\_147193) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GLIS1 (NM_147193) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GLIS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_147193
ORF Size:	1860 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211202).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_147193.1</a>
RefSeq Size:	2816 bp
RefSeq ORF:	1863 bp
Locus ID:	148979
UniProt ID:	<a href="#">Q8NBF1</a>
Cytogenetics:	1p32.3
MW:	66 kDa
Gene Summary:	GLIS1 is a GLI (MIM 165220)-related Kruppel-like zinc finger protein that functions as an activator and repressor of transcription (Kim et al., 2002 [PubMed 12042312]).[supplied by OMIM, Mar 2008]



[View online »](#)