

Product datasheet for **RC211213L2V**

EXOC3L2 (NM_138568) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EXOC3L2 (NM_138568) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EXOC3L2
Synonyms:	XTP7
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_138568
ORF Size:	1227 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211213).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138568.2
RefSeq Size:	1719 bp
RefSeq ORF:	1230 bp
Locus ID:	90332
UniProt ID:	Q2M3D2
Cytogenetics:	19q13.32
MW:	45.7 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is upregulated by vascular endothelial growth factor A and interacts with exocyst complex component 4. The encoded protein may be part of an exocyst complex that plays a role in cell membrane dynamics. Mutations in this gene may be associated with Alzheimer's disease. [provided by RefSeq, May 2017]