

Product datasheet for **RC211226L4V**

MED9 (NM_018019) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MED9 (NM_018019) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MED9
Synonyms:	MED25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018019
ORF Size:	438 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211226).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018019.2
RefSeq Size:	2222 bp
RefSeq ORF:	441 bp
Locus ID:	55090
UniProt ID:	Q9NWA0
Cytogenetics:	17p11.2
MW:	16.4 kDa



[View online »](#)

Gene Summary:

The multiprotein Mediator complex is a coactivator required for activation of RNA polymerase II transcription by DNA bound transcription factors. The protein encoded by this gene is thought to be a subunit of the Mediator complex. This gene is located within the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]