

Product datasheet for **RC211262L1V**

C20orf19 (KIZ) (NM_018474) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C20orf19 (KIZ) (NM_018474) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C20orf19
Synonyms:	C20orf19; HT013; Kizuna; NCRNA00153; PLK1S1; RP69
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_018474
ORF Size:	2019 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211262).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018474.3
RefSeq Size:	2329 bp
RefSeq ORF:	2022 bp
Locus ID:	55857
UniProt ID:	Q2M2Z5
Cytogenetics:	20p11.23
MW:	75.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene localizes to centrosomes, strengthening and stabilizing the pericentriolar region prior to spindle formation. The encoded protein usually remains with the mother centrosome after centrosomal duplication. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2013]