

Product datasheet for **RC211267A1V**

Human ARFIP1 (NM_001025593) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ARFIP1 (NM_001025593) AAV Particle
Tag:	Myc-DDK
Symbol:	ARFIP1
Synonyms:	HSU52521
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211267 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTCAAGAATCTCCCAAAATTCAGCAGCAGAAATCCAGTGACTAGTAATGGAGAAGTTGATGACT
CTCGTGAACATAGCTTTAATAGGGATTGAAGCATTATTACCATCTGGACTTGGTCTCTCAGAAACCCA
AATTACATCTCATGGCTTTGACAATACCAAGAGGGTGTATTGAAGCAGGAGCATTTCAAGGTGCCAG
AGAACACAGACAAAAAGTGACCAGTTATTCTAGCAGATGAAATTAATAATCCTGCAATGGAAAAAGTTAG
AACTTGTAGAAAAAGGAGTCTAAACACCTATAAGTGTACTCGACAGATTATCTCTGAGAAGCTAGGCCG
TGGCTCAAGAACTGTGGACCTTGAAGCTGAGATTGATATATTAAAGGGATAACAAGAAAAAATAT
GAAAAATTTTTAAAGCTGGCTCAAACTTGTGACCCAGCTTTTCCAGATGGTACATACCAAGGCAAC
TTGGAGATGCATTTGCTGACCTGAGTTTGAAGTCACTAGAACTTCATGAAGAATTTGGCTATAATGCCGA
TACCCAGAACTGCTGGCTAAAAATGGAGAGACTCTTCTTGGGGCCATTAATTTTTTTCATTGCTAGTGTG
AACACTTTGGTGAATAAAACATTGAAGATACATTAATGACTGTGAAACAGTATGAAAGTGCCAGGATTG
AATATGATGCATATCGCACTGATTTGGAAGAACTGAATCTTGGACCACGTGACGCAAACTCTGCCAA
GATTGAGCAGTCACAGCATCTTCCAAGCACATAAGGAAAAATATGATAAAATGCGCAATGATGTTCT
GTCAAATTGAAATTTCTAGAAGAAAATAAGGTTAAAGTATTGCACAATCAGCTGGTCCTTTCCACAATG
CCATTGCCGCTTACTTTGCTGGGAATCAGAAGCAGCTTGAACAGACACTTAACAGTTCCATATCAAAT
GAAACCCCTGGAGTGGATGCCCATCTTGGCTTGAAGAACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211267 protein sequence Red=Cloning site Green=Tags(s) MAQESPKNSAAEIPVTSNGEVDDSRHSFNRLKHSLPSGLGLSETQITSHGFDNTKEGVIEAGAFQGGQ RTQTKSGPVILADEIKNPAMEKLELVRKWSLNTYKCTRQIISEKLGRGSRTVDLELEAQIDILRDNKKKY ENILKLAQTLSTQLFQMVHTQRQLGDAFADLSLKSLELHEEFGYNADTQKLLAKNGETLLGAINFFIASV NTLVNKTIEDTLMTVKQYESARIEYDAYRTDLEELNLGPRDANTLPKIEQSQHLFQAHKEKYDKMRNDVS VKLKFLEENKVKVLHNQLVLFHNAIAAYFAGNQKQLEQTLKQFHIKLTGVDAPSWLEEQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001025593
ORF Size:	1023 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001025593.1, NP_001020764.1</u>
RefSeq Size:	3031 bp
RefSeq ORF:	1026 bp
Locus ID:	27236
UniProt ID:	<u>P53367</u>
Cytogenetics:	4q31.3
MW:	38.6 kDa