

Product datasheet for **RC211270A1V**

Human BAFF Receptor (TNFRSF13C) (NM_052945) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human BAFF Receptor (TNFRSF13C) (NM_052945) AAV Particle
Tag:	Myc-DDK
Symbol:	BAFF Receptor
Synonyms:	BAFF-R; BAFFR; BROMIX; CD268; CVID4; prolixin
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211270 representing NM_052945 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGCGAGGGCCCCGAGCCTGCGGGGACGGGACGCGCCAGCCCCACGCCCTGCGTCCCGGCCGAGT
GCTTCGACCTGCTGGTCCGCACTGCGTGGCCTGCGGGCTCCTGCGCACGCCGCGGCCGAAACCGGCCG
GGCCAGCAGCCCTGCGCCAGGACGGCGCTGCAGCCGAGGAGTCGGTGGGCGCGGGGGCCGGCAGGGCG
GCGCTGCCCTGCCGGGCTGCTCTTTGGCGCCCCGCGCTGCTGGGCTGGCACTGGTCTGCGCTGG
TCCTGGTGGTCTGGTGGTGGAGGCGGACAGCGGCGGCTTCGCGGCGGTCCTCCGAGAGGCCCC
CGACGGAGACAAGGACGCCCCAGAGCCCTGGACAAGGTATCATCTCTGTCTCCGGAATCTCTGATGCC
ACAGCTCCTGCCTGGCTCCTCCTGGGGAAGACCCAGGAACACCCACCTGGCCACAGTGTCCCTGTGC
CAGCCACAGAGCTGGGCTCCACTGAAGTGGTACCACCAAGACGGCCGGCCCTGAGCAACAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC211270 representing NM_052945 Red=Cloning site Green=Tags(s)
-------------------	--

MRRGPRSLRGRDAPAPTPCPVAECFDLLVRHCVACGLLRTPRKPAGASSPAPRTALQPQESVGAGAGEA
ALPLPGLLFGAPALLGLALVLALVLVGLVSWRRRQRRLRGASSAEAPDGDKDAPEPLDKVILSPGISDA
TAPAWPPPGEDPGTTPPGHVPVPATELGSTELVTTKTAGPEQQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species:	Human
Serotype:	AAV-2



[View online »](#)

ACCN:	NM_052945
ORF Size:	552 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_052945.2</u>
RefSeq Size:	898 bp
RefSeq ORF:	555 bp
Locus ID:	115650
UniProt ID:	<u>Q96R13</u>
Cytogenetics:	22q13.2
MW:	18.7 kDa