

## Product datasheet for **RC211275L4V**

### SLURP1 (NM\_020427) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLURP1 (NM_020427) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLURP1
Synonyms:	ANUP; ARS; ArsB; LY6-MT; LY6LS; MDM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020427
ORF Size:	309 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211275).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020427.2</a>
RefSeq Size:	537 bp
RefSeq ORF:	312 bp
Locus ID:	57152
UniProt ID:	<a href="#">P55000</a>
Cytogenetics:	8q24.3
Protein Families:	Druggable Genome, Secreted Protein
MW:	11.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the Ly6/uPAR family but lacks a GPI-anchoring signal sequence. It is thought that this secreted protein contains antitumor activity. Mutations in this gene have been associated with Mal de Meleda, a rare autosomal recessive skin disorder. This gene maps to the same chromosomal region as several members of the Ly6/uPAR family of glycoprotein receptors. [provided by RefSeq, Jul 2008]