

Product datasheet for **RC211295A1V**

Human ODF4 (NM_153007) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ODF4 (NM_153007) AAV Particle
Tag:	Myc-DDK
Symbol:	ODF4
Synonyms:	CT134; CT136; OPPO1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211295 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATGCAGAGTACTCTGGGAATGAGTTCCCCAGGTCAGAAGGAGAAAGAGACCAACATCAGAGACCTG
GAAAGGAAAGGAAGAGTGGGGAGGCAGGACGGGGCACAGGTGAGCTGGGACAAGATGGGAGACTGCTGTC
CTCCACCCTCTCCCTCAGTAGTAACAGGTCCTTGGGCCAGCGCCAGAAGTCTCCGCTGCCCTTTCAATGG
AGAATCACACACAGCTTCCGCTGGATGGCCAGGTGTTGGCCTCTGAGCTCAGCCTGGTTGCCTTTATCC
TACTATTGGTCATGGCCTTCTCCAAGAAATGGCTGGACCTCTCTAGGAGCCTCTTCTACCAGCGCTGGCC
CGTGGATGTGAGCAACAGAATCCACACATCAGCCACGTTATGTCCATGGGGCTCCTGCACCTTTTGCAAA
TCCAGGAGCTGTTCTGACTTAGAGAATGGGAAAGTCACCTTCATCTTCTCCACCCTCATGCTATTCCCCA
TTAACATCTGGATCTTCGAGTTGGAAAGGAATGTATCCATCCCCATAGGCTGGAGCTATTTATTGGTTG
GCTGGTGCTTATCCTATACTTACCTGCGCGATCCTTTGCTACTTCAACCATAAAAGTTTCTGGAGTCTG
ATTCTGAGCCACCCAGTGGTGGCGTGTCTGCAGCAGCAGTTTCGGCTCAGTAGAAGAATCTCCAAGGG
CACAGACGATCACAGACCCCCATCACCAGGAGGGAGTCCTGGATCCTGAGCAGAAGGATACACATGT
G

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211295 protein sequence Red=Cloning site Green=Tags(s) MDAEYSGNEFPRSEGERDQHQRPGKERKSGEAGRGTGELGQDGRLLSSTLSLSSNRSLGQRQNSPLPFQW RITHSFRWMAQVLASELSLVAFILLLVMAFSKKWLDLSRSLFYQRWPDVSNRIHTSAHVMSMGLLHFCK SRSCSDLENGKVTFIFSTLMLFPINIWIFELERNVSIPIGWSYFIGWLVLILYFTCAILCYFNHKSFWSL ILSHPSGAVSCSSSFGSVEESPRATITDTPITQEGVLDPEQKQDTHV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_153007
ORF Size:	771 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_153007.3</u>
RefSeq Size:	1149 bp
RefSeq ORF:	774 bp
Locus ID:	146852
UniProt ID:	<u>Q2M2E3</u>
Cytogenetics:	17p13.1
MW:	29.2 kDa